

Product Datasheet

Mouse Galectin 8 (GAL8) ELISA Kit (orb1207600)

Description	Mouse Galectin 8 (GAL8) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9JL15
Sensitivity	0.055ng/mL
Entrez	56048
Expiration Date	Please enquire.