

Product Datasheet

Mouse Fibrinogen Like Protein 2 (FGL2) ELISA Kit (orb1207613)

Description	Mouse Fibrinogen Like Protein 2 (FGL2) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12804
Sensitivity	5.4pg/mL
Entrez	14190
Expiration Date	Please enquire.