

Product Datasheet

Mouse Fatty Acid Binding Protein 6, Ileal (FABP6) ELISA Kit (orb1207616)

Description	Mouse Fatty Acid Binding Protein 6, Ileal (FABP6) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P51162
Sensitivity	23pg/mL
Entrez	16204
Expiration Date	Please enquire.