

Product Datasheet

Mouse Dihydrolipoyl Transacetylase (DLAT) ELISA Kit (orb1207628)

Description	Mouse Dihydrolipoyl Transacetylase (DLAT) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8BMF4
Sensitivity	0.048ng/mL
Entrez	235339
Expiration Date	Please enquire.