

Product Datasheet

Mouse Defensin Beta 2 (DEFb2) ELISA Kit (orb1207630)

Description	Mouse Defensin Beta 2 (DEFb2) ELISA Kit
Reactivity	Mouse
Range	3.12-200pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P82020
Sensitivity	1.31pg/mL
Entrez	13215
Expiration Date	Please enquire.