

Product Datasheet

Mouse Cyclophilin A (CYPA) ELISA Kit (orb1207633)

Description	Mouse Cyclophilin A (CYPA) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17742
Sensitivity	0.065ng/mL
Entrez	268373
Expiration Date	Please enquire.