

Product Datasheet

Mouse Cytochrome P450 2E1 (CYP2E1) ELISA Kit (orb1207634)

Description	Mouse Cytochrome P450 2E1 (CYP2E1) ELISA Kit
Reactivity	Mouse
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q05421
Sensitivity	1.29ng/mL
Entrez	13106
Expiration Date	Please enquire.