

Product Datasheet

Mouse Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit (orb1207636)

Description	Mouse Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P70658
Sensitivity	27.2pg/mL
Entrez	12767
Expiration Date	Please enquire.