

Product Datasheet

Mouse Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit (orb1207639)

Description	Mouse Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	0.046ng/mL
Expiration Date	Please enquire.