

Product Datasheet

Mouse Creatine Kinase, Muscle (CKM) ELISA Kit (orb1207657)

Description	Mouse Creatine Kinase, Muscle (CKM) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07310
Sensitivity	0.55ng/mL
Entrez	12715
Expiration Date	Please enquire.