

Product Datasheet

Mouse CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit (orb1207659)

Description	Mouse CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P28033
Sensitivity	33pg/mL
Entrez	12608
Expiration Date	Please enquire.