

Product Datasheet

Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit (orb1207663)

Description	Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P32972
Sensitivity	0.032ng/mL
Entrez	21949
Expiration Date	Please enquire.