

Product Datasheet

Mouse Carbonic Anhydrase I (CA1) ELISA Kit (orb1207674)

Description	Mouse Carbonic Anhydrase I (CA1) ELISA Kit
Reactivity	Mouse
Range	0.781-50IU/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13634
Sensitivity	0.3IU/mL
Entrez	12346
Expiration Date	Please enquire.