

Product Datasheet

Mouse Aquaporin 4 (AQP4) ELISA Kit (orb1207691)

Description	Mouse Aquaporin 4 (AQP4) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P55088
Sensitivity	23.7pg/mL
Entrez	11829
Expiration Date	Please enquire.