

Product Datasheet

Mouse Aquaporin 2, Collecting Duct (AQP2) ELISA Kit (orb1207693)

Description	Mouse Aquaporin 2, Collecting Duct (AQP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P56402
Sensitivity	0.059ng/mL
Entrez	11827
Expiration Date	Please enquire.