

Product Datasheet

Mouse Apolipoprotein E (APOE) ELISA Kit (orb1207696)

Description	Mouse Apolipoprotein E (APOE) ELISA Kit
Reactivity	Mouse
Range	7.81-500ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08226
Sensitivity	3.16ng/mL
Entrez	11816
Expiration Date	Please enquire.