

Product Datasheet

Mouse Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit (orb1207706)

Description	Mouse Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q3TIW9
Sensitivity	0.051ng/mL
Entrez	238011
Expiration Date	Please enquire.