

Product Datasheet

Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit (orb1207717)

Description	Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q7TMS5
Sensitivity	0.055ng/mL
Entrez	26357
Expiration Date	Please enquire.