

Product Datasheet

Human Proteasome 26S Subunit, Non ATPase 4 (PSMD4) ELISA Kit (orb1207735)

Description	Human Proteasome 26S Subunit, Non ATPase 4 (PSMD4) ELISA Kit
Reactivity	Human
Range	78-5000pg/ml
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P55036
Sensitivity	29.8pg/mL
Entrez	5710
Expiration Date	Please enquire.