

Product Datasheet

Human X-Ray Repair Cross Complementing 6 (XRCC6) ELISA Kit (orb1207738)

Description	Human X-Ray Repair Cross Complementing 6 (XRCC6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P12956
Sensitivity	0.112ng/mL
Entrez	2547
Expiration Date	Please enquire.