

Product Datasheet

Human Xeroderma Pigmentosum, Complementation Group C (XPC) ELISA Kit (orb1207739)

Description	Human Xeroderma Pigmentosum, Complementation Group C (XPC) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q01831
Sensitivity	0.115ng/mL
Entrez	7508
Expiration Date	Please enquire.