

Product Datasheet

Human WAP Four Disulfide Core Domain Protein 2 (WFDC2) ELISA Kit (orb1207740)

Description	Human WAP Four Disulfide Core Domain Protein 2 (WFDC2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, saliva, tissue homogenates or other biological fluids.
Uniprot ID	Q14508
Sensitivity	0.068ng/mL
Entrez	10406
Expiration Date	Please enquire.