

Product Datasheet

Human Von Willebrand Factor A Domain Containing Protein 3A (vWA3A) ELISA Kit (orb1207742)

Description	Human Von Willebrand Factor A Domain Containing Protein 3A (vWA3A) ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	plasma.
Uniprot ID	A6NCI4
Sensitivity	0.15ng/mL
Entrez	146177
Expiration Date	Please enquire.