

Product Datasheet

Human Vasoactive Intestinal Peptide (VIP) ELISA Kit (orb1207743)

Description	Human Vasoactive Intestinal Peptide (VIP) ELISA Kit
Reactivity	Human
Range	7.81-500pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01282
Sensitivity	2.34pg/mL
Entrez	7432
Expiration Date	Please enquire.