

Product Datasheet

Human Toll Like Receptor 6 (TLR6) ELISA Kit (orb1207768)

Description	Human Toll Like Receptor 6 (TLR6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Y2C9
Sensitivity	0.112ng/mL
Entrez	10333
Expiration Date	Please enquire.