

Product Datasheet

Human TGF Beta Inducible Early Response Gene 1 (TIEG1) ELISA Kit (orb1207774)

Description	Human TGF Beta Inducible Early Response Gene 1 (TIEG1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13118
Sensitivity	0.122ng/mL
Entrez	7071
Expiration Date	Please enquire.