

Product Datasheet

Human Alpha-Fodrin (SPTAN1) ELISA Kit (orb1207788)

Description	Human Alpha-Fodrin (SPTAN1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q13813
Sensitivity	0.064ng/mL
Entrez	6709
Expiration Date	Please enquire.