

Product Datasheet

Human Surfactant Associated Protein A (SPA) ELISA Kit (orb1207789)

Description	Human Surfactant Associated Protein A (SPA) ELISA Kit
Reactivity	Human
Range	1.87-120ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
Uniprot ID	Q8IWL2
Sensitivity	0.69ng/mL
Entrez	653509
Expiration Date	Please enquire.