

Product Datasheet

Human Sex Determining Region Y Box Protein 1 (SOX1) ELISA Kit (orb1207790)

Description	Human Sex Determining Region Y Box Protein 1 (SOX1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	000570
Sensitivity	0.059ng/mL
Entrez	6656
Expiration Date	Please enquire.