

Product Datasheet

Human Superoxide Dismutase 3, Extracellular (SOD3) ELISA Kit (orb1207793)

Description	Human Superoxide Dismutase 3, Extracellular (SOD3) ELISA Kit
Reactivity	Human
Range	6.25-400ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	P08294
Sensitivity	2.23ng/mL
Entrez	6649
Expiration Date	Please enquire.