

Product Datasheet

Human Protein Kinase, AMP Activated Alpha 1 (PRKAa1) ELISA Kit (orb1207831)

Description	Human Protein Kinase, AMP Activated Alpha 1 (PRKAa1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13131
Sensitivity	0.062ng/mL
Entrez	5562
Expiration Date	Please enquire.