

Product Datasheet

Human Peroxisome Proliferator Activated Receptor Alpha (PPARα) ELISA Kit (orb1207834)

Description	Human Peroxisome Proliferator Activated Receptor Alpha (PPARα) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q07869
Sensitivity	0.052ng/mL
Entrez	5465
Expiration Date	Please enquire.