

Product Datasheet

Human Programmed Cell Death Protein 1 Ligand 2 (PDCD1LG2) ELISA Kit (orb1207866)

Description	Human Programmed Cell Death Protein 1 Ligand 2 (PDCD1LG2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BQ51
Sensitivity	0.063ng/mL
Entrez	80380
Expiration Date	Please enquire.