

Product Datasheet

Human Nitric Oxide Synthase 1, Neuronal (NOS1) ELISA Kit (orb1207898)

Description	Human Nitric Oxide Synthase 1, Neuronal (NOS1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P29475
Sensitivity	0.053ng/mL
Entrez	4842
Expiration Date	Please enquire.