

Product Datasheet

Human Mannose Receptor C Type 2 (MRC2) ELISA Kit (orb1207916)

Description	Human Mannose Receptor C Type 2 (MRC2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9UBG0
Sensitivity	0.117ng/mL
Entrez	9902
Expiration Date	Please enquire.