

Product Datasheet

Human Intercellular Adhesion Molecule 3 (ICAM3) ELISA Kit (orb1207984)

Description	Human Intercellular Adhesion Molecule 3 (ICAM3) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P32942
Sensitivity	0.4ng/mL
Entrez	3385
Expiration Date	Please enquire.