

Product Datasheet

Human Heme Oxygenase 2, Decycling (HO2) ELISA Kit (orb1207998)

Description	Human Heme Oxygenase 2, Decycling (HO2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P30519
Sensitivity	0.066ng/mL
Entrez	3163
Expiration Date	Please enquire.