

Product Datasheet

Human Glutathione S Transferase Alpha 3 (GSTa3) ELISA Kit (orb1208017)

Description	Human Glutathione S Transferase Alpha 3 (GSTa3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q16772
Sensitivity	0.053ng/mL
Entrez	2940
Expiration Date	Please enquire.