

Product Datasheet

Human Gonadotropin Releasing Hormone (GnRH) ELISA Kit (orb1208024)

Description	Human Gonadotropin Releasing Hormone (GnRH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01148
Sensitivity	0.064ng/mL
Entrez	2796
Expiration Date	Please enquire.