

Product Datasheet

Human Glucosidase Alpha, Acid (GaA) ELISA Kit (orb1208045)

Description	Human Glucosidase Alpha, Acid (GaA) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P10253
Sensitivity	0.064ng/mL
Entrez	2548
Expiration Date	Please enquire.