

Product Datasheet

Human Delta-Sleep Inducing Peptide (dSIP) ELISA Kit (orb1208071)

Description	Human Delta-Sleep Inducing Peptide (dSIP) ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q99576
Sensitivity	47pg/mL
Entrez	1831
Expiration Date	Please enquire.