

Product Datasheet

Human Defensin Beta 2 (DEFb2) ELISA Kit (orb1208080)

Description	Human Defensin Beta 2 (DEFb2) ELISA Kit
Reactivity	Human
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	015263
Sensitivity	13.55pg/mL
Entrez	1673
Expiration Date	Please enquire.