

Product Datasheet

Human Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit (orb1208089)

Description	Human Chemokine C-X-C-Motif Receptor 3 (CXCR3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P49682
Sensitivity	0.061ng/mL
Entrez	2833
Expiration Date	Please enquire.