

Product Datasheet

Human Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit (orb1208093)

Description	Human Cross Linked C-Telopeptide Of Type II Collagen (CTXII) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	0.045ng/mL
Expiration Date	Please enquire.