

Product Datasheet

Human Chorionic Somatomammotropin Hormone 1 (CSH1) ELISA Kit (orb1208100)

Description	Human Chorionic Somatomammotropin Hormone 1 (CSH1) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	PODML2
Sensitivity	1.12ng/mL
Entrez	1442
Expiration Date	Please enquire.