

Product Datasheet

Human Creatine Kinase, Mitochondrial 1A (CKMT1A) ELISA Kit (orb1208110)

Description	Human Creatine Kinase, Mitochondrial 1A (CKMT1A) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P12532
Sensitivity	0.57ng/mL
Entrez	548596
Expiration Date	Please enquire.