

Product Datasheet

Human Creatine Kinase, Muscle (CKM) ELISA Kit (orb1208111)

Description	Human Creatine Kinase, Muscle (CKM) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06732
Sensitivity	0.63ng/mL
Entrez	1158
Expiration Date	Please enquire.