

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Alpha 1 (CHRNA1) ELISA Kit (orb1208112)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 1 (CHRNA1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P02708
Sensitivity	0.28ng/mL
Entrez	1134
Expiration Date	Please enquire.