

Product Datasheet

Human CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit (orb1208119)

Description	Human CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P17676
Sensitivity	0.061ng/mL
Entrez	1051
Expiration Date	Please enquire.