

Product Datasheet

Human CCAAT/Enhancer Binding Protein Alpha (CEBPa) ELISA Kit (orb1208120)

Description	Human CCAAT/Enhancer Binding Protein Alpha (CEBPa) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P49715
Sensitivity	0.055ng/mL
Entrez	1050
Expiration Date	Please enquire.